


EF Series



**130°
Bending**



**1280*800
Resolution**



Specifications	EF20	EF28	EF30	EF38	EF45	EF52
Outer diameter of insertion tube	Φ2.0mm	Φ2.8mm	Φ3.0mm	Φ3.8mm	Φ4.5mm	Φ5.2mm
Inner diameter of instrument channel	/	/	Φ1.2mm	/	Φ2.2mm	Φ2.6mm
Field of view	120°					
Depth of view field	3mm-200mm					
Range of tip Bending	(U)130° (D)130°					
Working length	360mm					

Specifications	VS100
Imaging resolution	1280x800
Video output	HDMI
Data storage	USB Disk, SD card(Built-in 8GB)
Inner/ Video format	JPG/ BMP/ MP4
Power	DC3.7V Li-ion Battery; 100~240V; 50/60Hz

-  +86 400-0756-365 (Service)
-  +86 (0) 756 6910150
-  301 Building 4, No. 30 Feipeng Road, GMTCM Park,
Hengqin New Area, Zhuhai City, Guangdong Province, China
-  international@vision-medt.com



LINKEDIN
Vision-medt



TWITTER
Vision-medt

www.vision-medt.com